

DATA SHEET

Product overview

Name CELT-300

Short description Potent and selective hA_{2A} Adenosine receptor fluorescent ligand

Biological description It shows full selectivity for A_{2A} over A₁, A_{2B} and A₃ (only in the A_{2A}

receptor it is possible to measure a K_i whose value is 8.35 nM in

radioligand binding assay).

Biological action Selective orthosteric ligand

Quantity 10 μg

Purity > 97%

Properties

Molecular Weight 1320.53 (TFA Salt)

Source Synthetic

Appearance Dark blue solid

Formulation Solid powder

Excitation 646 nm

Emission 662 nm

Pharmacological validation The efficacy and potency of CELT-300 as a selective fluorescent hA_{2A}

adenosine receptor ligand was confirmed by a radioligand binding

assay.

Storing and Using product

Storage instructions -20 °C (protect from light)

Solubility overview Soluble in DMSO

Stock solution Add 76 µL of assay buffer containing 1% of DMSO to obtain a 100 µM

stock solution.

Handling After thawing individual aliquots for use, we recommend briefly

sonicating the sample to ensure it is fully dissolved and the solution is

homogeneous. We do not recommend using the product after subjecting it to repetitive freeze-thaw cycles.

Shipping conditions The product, as a solid, is stable at ambient temperature for periods

of up to a few days and does not require shipping on ice/dry ice.

Important This product is for RESEARCH USE ONLY and is not intended for

therapeutic or diagnostic use. Not for human or veterinary use.