



**Celtarys**  
RESEARCH



# **CELT-150 VHL E3 Ligase** **TR-FRET protocol**

# Protocol for E3 ubiquitin ligase VHL affinity binding assay using **TR-FRET**

CELT-150, a potent VHL E3 ligase fluorescent ligand, was used as fluorescent probe in a TR-FRET assay to study the affinity of a known inhibitor of E3 ubiquitin ligase VHL, VH298.

## Equipment, Materials and Reagents:

- Plate Reader with a TR-FRET optic module capable of emitting in the 390nm range and detecting in the 650nm range.
- TR-FRET compatible plates such as *CulturPlate-384*.
- Centrifuge with adaptor compatible with 384 plates.
- Recombinant His Tagged Human VBC complex.
- Anti-His Europium labeled Antibody.
- CELT-150, VHL E3 ligase fluorescent ligand

### Assay Buffer:

- 20mM HEPES pH 7.5
- 150mM NaCl
- 1mM TCEP
- 0.01% TWEEN20

### Master Mix:

In order to streamline the process and run as many compounds as possible, preparation of a master mix is recommended. This mixture should have 2 times the concentration of the final concentration in the well:

- 400nM CELT-150.
- 100nM His-VBC complex.
- 1nM AntiHis-EuAb.
- Assay Buffer.

## Assay Procedure

The assay was performed using 384 Microplates with 15uL as a final assay volume.

The assay was performed using Europium (Anti-His-EuAb) as donor and His-VBC complex as protein of interest to be degraded.

A 2X compound concentration template plate is recommended diluting the compounds in the assay buffer. The final DMSO concentration in the wells should be 2%.

7.5uL of each 2X template well are deposited in their corresponding wells, and once all compounds are ready the plate is centrifuged.

Afterwards, 7.5uL of the Master Mix are added to each well, totaling 15uL in each well. The plate is centrifuged again.

It is recommended to cover the plate from the light while incubating. Measures were taken every half an hour for 3hours, as the signal stabilizes around the 90minute mark.


## Assay Performance

·Z´ prime of 0.87.

·Different percentages of DMSO were tested, from 0.1% to 2% (final concentration in the wells). Up to 2% is acceptable although the assay window is reduced with the increasing percentages of DMSO.

## Results

The figure below shows how increasing concentrations of VH298 compete with CELT-150 ligand and thereby prevent TR-FRET from occurring.



VH298 displayed the expected potency in good correlation with the literature (Frost et al., 2016, Potent and selective chemical probe of hypoxic signalling downstream of HIF- $\alpha$  hydroxylation via VHL inhibition. Nature Communications volume 7, Article number: 13312).