



## DATA SHEET

### Product overview

Name	CELT-300
Short description	Potent and selective hA <sub>2A</sub> Adenosine receptor fluorescent ligand
Biological description	It shows full selectivity for A <sub>2A</sub> over A <sub>1</sub> , A <sub>2B</sub> and A <sub>3</sub> (only in the A <sub>2A</sub> receptor it is possible to measure a K <sub>i</sub> whose value is 8.35 nM in radioligand binding assay).
Biological action	Selective orthosteric ligand
Quantity	50 µg
Purity	> 97%

---

### Properties

Molecular Weight	1334.42
Source	Synthetic
Appearance	Dark blue solid
Formulation	Solid powder
Excitation	646 nm
Emission	662 nm

---

### Applications

Pharmacological validation	The efficacy and potency of CELT-300 as a selective fluorescent hA <sub>2A</sub> adenosine receptor ligand. was confirmed by a radioligand binding assay.
----------------------------	---

---

### Storing and Using product

Storage instructions	-20 °C (protect from light)
Solubility overview	Soluble in DMSO
Stock solution	Add 37 µL of DMSO to obtain a 1mM solution

Handling	After thawing individual aliquots for use, we recommend briefly sonicating the sample to ensure it is fully dissolved and the solution is homogeneous. We do not recommend using the product after subjecting it to repetitive freeze-thaw cycles.
Shipping conditions	The product, as a solid, is stable at ambient temperature for periods of up to a few days and does not require shipping on ice/dry ice.
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.