



## DATA SHEET

### Product overview

Name	CELT-429
Short description	Potent and selective hD <sub>3</sub> Dopamine receptors fluorescent ligand
Biological description	It shows full selectivity for D <sub>3</sub> over D <sub>2</sub> and D <sub>4</sub> (only in the D <sub>3</sub> receptor it is possible to measure a K <sub>i</sub> whose value is 75.4 nM) in a radioligand binding assay.
Biological action	Selective orthosteric ligand
Quantity	50 µg
Purity	> 97%

---

### Properties

Molecular Weight	1277.13
Source	Synthetic
Appearance	Purple solid
Formulation	Lyophilized solid
Excitation	589 nm
Emission	616 nm
Pharmacological validation	The efficacy and potency of CELT-429 as a selective hD <sub>3</sub> fluorescent ligand was confirmed by a radioligand binding assay.

---

### Storing and Using product

Storage instructions	-20 °C (protect from light)
Solubility overview	Soluble in DMSO.
Stock solution	Add 39 µL of DMSO to obtain a 1mM solution
Handling	After thawing individual aliquots for use, we recommend briefly sonicating the sample to ensure it is fully dissolved and the solution is homogeneous. We do not recommend using the product after subjecting it to repetitive freeze-thaw cycles.

Shipping conditions

The product, as a solid, is stable at ambient temperature for periods of up to a few days and does not require shipping on ice/dry ice.

Important

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.